



ABIF workshop

Introduction to writing macros in Fiji

with Joel Ryan
McGill University
21.10.2021

Introduction to writing macros in Fiji

- Goal:
 - Learn a general workflow to write macros in Fiji in order to help:
 - Process images
 - Extract data from images
- You will need to download:
 - Fiji
 - fiji.sc
 - Sample image data:
 - https://github.com/ABIF-McGill/Intro_to_Fiji_Macros



Why write macros/scripts?

- Streamline processing and data extraction
- Help develop more robust and reliable assays

--> Relieve the bottleneck of tedious analysis

Open image, adjust contrast, duplicate


Apply convolutional filters

Calculate and apply a threshold

Generate ROIs, select target image, apply ROIs to other channel(s)

Measure intensity data in each ROI

Repeat this process on
10-30 images x 3 replicates x n genotypes/conditions...



Today's workshop

- Walk through writing three macros
 1. Make a composite image
 2. Apply filters, threshold and measure intensities
 3. Manual analysis template / quality control

Today's workshop

- Walk through writing three macros
 1. Make a composite image
 2. Apply filters, threshold and measure intensities
 3. Manual analysis template / quality control
- > IJ1 macro language
 - Targeted for image processing/analysis
 - Written commands underlying point-and-click ImageJ/Fiji
 - Great intro to programming for bio people

Today: we are a big group!

- There are lots of participants today!
- To keep things moving:
 - Please post questions in the Zoom chat
 - We'll try our best to answer as many questions as we can
- If things aren't working:
 - Take note of the error message and line
 - Try copy-pasting the script from the github repository...

Challenges for you beyond this workshop


- Figure out the processing and analysis needed for your imaging experiments:
 - Conceptually, what is the question?
 - -> practically:
 - what quantitative parameters do you need to extract from your images
 - which functions do you need to run to get that data
 - do you need to retool your image acquisition settings...?
- ... then find good ways to plot and present your data!

Exercise 1 – raw vs composite RGB images



demo_DAPI_Phalloidin_Mitotracker_001.tif

Raw image, drag-and-drop in powerpoint.



demo_DAPI_Phalloidin_Mitotracker_001.tif (RGB).tif

Composite RGB image generated with Exercise 1 macro,
drag-and-drop in powerpoint.

for loop

```
for (i = 0; i < numROI; i++) {  
    roiManager("Select", i);  
    Stack.setChannel(2);  
    run("Measure");  
}
```

for loop

Start value

End value

Increment

```
for (i = 0; i < numROI; i++) {  
  
    roiManager("Select", i);  
    Stack.setChannel(2);  
    run("Measure");  
  
}
```

for loop

Start value

End value

Increment

first, “build a sequence of integers going from the Start Value, to the End value, with an increment of 1”

```
for (i = 0; i < numROI; i++) {  
  
    roiManager("Select", i);  
    Stack.setChannel(2);  
    run("Measure");  
  
}
```

for loop

Start value

End value

Increment

```
for (i = 0; i < numROI; i++) {  
  
    roiManager("Select", i);  
    Stack.setChannel(2);  
    run("Measure");  
  
}
```

first, “build a sequence of integers going from the Start Value, to the End value, with an increment of 1”

---in the case of numROI being 5:

0, 1, 2, 3, 4

for loop

Start value

End value

Increment

```
for (i = 0; i < numROI; i++) {  
  
    roiManager("Select", i);  
    Stack.setChannel(2);  
    run("Measure");  
  
}
```

first, “build a sequence of integers going from the Start Value, to the End value, with an increment of 1”

---in the case of numROI being 5:

0, 1, 2, 3, 4,

“Run the { content of the loop } replacing i with the first integer of the sequence.

Then run the { content of the loop } replacing i with the second integer in the sequence... “ and so on...

for loop

Start value End value Increment

```
for (i = 0; i < numROI; i++) {  
  
    roiManager("Select", i);  
    Stack.setChannel(2);  
    run("Measure");  
  
}
```

```
i = 0  
roiManager("Select", i);  
Stack.setChannel(2);  
run("Measure");
```

```
i = 1  
roiManager("Select", i);  
Stack.setChannel(2);  
run("Measure");
```

```
i = 2  
roiManager("Select", i);  
Stack.setChannel(2);  
run("Measure");
```

```
i = 3  
roiManager("Select", i);  
Stack.setChannel(2);  
run("Measure");
```

...